[PACKAGING SPECIFICATION]

Variant	1 piece in a pouch	5 pieces in a pouch	5 pieces in a box	25 pieces in a package	50 pieces in a package
L-CAB-Pro/1B - 1P	1	5	25	50	
L-CAB-Pro/5B			25	50	
L-CAB-Pro/25B				50	
L-CAB-Pro/50B					50

[INTRODUCTION]

The novel coronaviruses belong to the *Beta* genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and headache are also common. COVID-19 is highly contagious and can cause viral infection of the respiratory tract. It is an infectious disease, easily spreading through coughing and sneezing droplets containing viruses. Type A viruses are typically more prevalent than type B viruses and are associated with most severe COVID-19 epidemics. Type A viruses are usually detected in hospital diagnosis or isolation. Type B has become more important due to the availability of effective antiviral therapy. Rapid diagnosis of influenza can reduce hospital stays, antimicrobial use and cost of hospital care.

[INTENDED USE]

LOMINA SARS-CoV-2 + Influenza A/B Antigen Test is an antigen test kit designed for simultaneous qualitative detection of the SARS-CoV-2 and/or Influenza A/B infection in human nasal swab samples in vitro. It is used as a supplementary detection indicator for suspected cases of SARS-CoV-2 causing disease COVID-19 and/or respiratory Flu infection A/B.

[TEST PRINCIPLE]

The Lomina SARS-CoV-2 + Influenza A/B Antigen Test is a qualitative immunoassay for the detection of SARS-CoV-2 and/or Flu A/ B antigen in nasal swab specimens. Anti-SARS-CoV-2 antibody is coated in the test region of the test card. The anti-Influenza A/B antibody is coated in the test A region (A) and test B region (B) respectively. During testing, antigen in the specimen reacts with anti-SARS-CoV-2 antibody-coated particles and with anti-Influenza A antibody-coated particles as well as with anti-Influenza B antibody-coated particles. The reaction pad is produced by immunological complex formation. The complex migrates along the membrane by capillary action to the test region. The complex then respectively reacts with the anti-SARS-CoV-2 antibody in the test region of SARS-CoV-2 and with anti-Influenza A antibody in the A region as well as with anti-Influenza B antibody in the B region. If the complex is positive, a colored line will appear in the test region of SARS-CoV-2. If the specimen contains Influenza A/B, a colored line will appear in the test A/B region. The specimen does not contain antigen of SARS-CoV-2, Influenza A/B, no colored line will appear in the test region, indicating a negative result. This serves as a control, a colored line will always appear in the control region (C), indicating that the result is effective.

[CONTENT OF PRESENTATION]

1. IVD test strip in a plastic cassette.
2. Humidity absorbing pad.
3. Stabilization fluid (Buffer).
4. Nasal Swab for sample collection **CE 0197**.
5. Driptip (100 µl).

6. Bottles containing 300 µl (+/-5%) of buffer that are being used as extraction tubes.
7. Instruction for Use.

[NOT INCLUDED IN THE PACKAGE]

Stopwatch/clock for measuring the duration of the test.
[PRECAUTIONS]

- The LOMINA SARS-CoV-2+Influenza A/B Antigen Test is intended for professional in vitro diagnostic use and should only be used for the qualitative detection of the presence of SARS-CoV-2 antigen and/or

- Read the Package Insert (IFU) prior to use. Instructions should be read and followed carefully.
- Do not use kit or components beyond the expiration date.
- The device contains material of animal origin and should be handled as a potential hazard.
- All devices are packed in foil pouches that exclude moisture during storage. Inspect each foil pouch before opening. Do not use devices that have holes in the foil or where the pouch has not been completely sealed. Erroneous result may occur if test reagents or components are inadvertently discarded.
- Do not use the Sample Diluent Buffer if it is discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
- All patient specimens should be handled and discarded as if they were biologically hazardous. Prior to testing, all samples must be thoroughly mixed with the reagent to ensure the correct solution for testing.
- Do not mix different lots (LOT) of tests and reagents!

[STORAGE CONDITIONS]

1. The test package shall be stored at 4°C to 30°C.
2. Kit contents are stable until the expiration dates marked on its outer packaging containers, its validity period is set to 24 months (based on stability study).

3. Prolonged exposure to heat and humidity will make the reagent useless.

[TRANSPORT CONDITIONS]

The test package may be temporarily transported at -10°C.

[SPECIMEN COLLECTION AND STORAGE]

Specimen collection:
• Acceptable specimens for testing with the LOMINA SARS-CoV-2+Influenza A/B Antigen Test include samples from nasal swabs.

• Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of sample with the interpretation of test results.

• Use freshly collected samples for best test results.

• To ensure optimal performance, use ONLY THOSE SWABS supplied in the kit.

• THE TEST IS DESIGNED TO DETECT AN ACTIVE VIRUS! Inactivated virus may not be detected / the test sensitivity will be affected!

• Do not test nasal swab specimens intended for RT-PCR, PCR or any other chemically treated.

• After sampling, the appropriate buffer provided in the kit shall be used as soon as possible. No other type of solution should be used to store the sample.

• Do not test preserved or otherwise prepared samples intended for PCR, ELISA and other methods!

[ANTERIOR NASAL SWAB PROCEDURE]

Use a nasal swab to take a sample. Insert the swab on the picture, carefully and slowly into the nostril. Push the swab until it meets resistance at the wall of the turbinates. Rotate the swab gently against the nasal wall a few times. If coughing, please wait a while and try again. After collection place the swab in the reagent bottle containing buffer.

If the swab body extends beyond the upper level of the reagent solution, squeeze it so that the top of the swab stem is just below the upper level of the solution and allow the sample to mix sufficiently with the buffer.

If the swab body extends beyond the upper level of the reagent solution, squeeze it so that the top of the swab stem is just below the upper level of the solution and allow the sample to mix sufficiently with the buffer.

When pushing/pulling the swab, keep it gentle and do not use force to avoid injury to the patient!

[SAMPLE QUALITY REQUIREMENT]

1. The sample should be processed by preservation solution (buffer) immediately after sample collection. It is recommended to detect directly without sample collection.

2. Reagents should not be opened until sample is applied.

3. Long-term storage of the samples is not recommended.

[OPERATING STEPS]

1. Before testing, read the operating instructions carefully, and restore the testing kit.

2. Take out the test card in the environment of temperature between 10°C and 30°C and humidity 30% to 50% (avoid strong convection or ventilation environment).

3. Collect the sample using SWAB as described above.

4. Place the SWAB with the swab into the tube with buffer and wash the tube well by rotating the SWAB at least 5 times, then remove the SWAB.


5. Close the tube by a dripping lid.

6. Tear open the foil bag, take out the test card, place on a horizontally leveled table and drip 3 drops (about 100 µl) of the treated sample solution into the sample hole of the test card.

7. Use the test card as soon as possible but no later than after 15 min. after

[INTERPRETATION OF TEST RESULTS]

The test should be judged in a combination with clinical symptoms and other detecting indicators. Result is interpreted by presence of a visible strip at the „COV-2, Flu A or B“ area at the test card body (see next picture), the visual inspection results are following:

Opening the pouch. Humidity may damage the test sensitivity. Apply the seal and wait for 15 minutes without moving the card. The test is to be invalid after more than 25 minutes after dripping the analyte (sample) in to detection card.

8. Insert to a depth of 2-3 cm

9. Apply the seal and wait for 15 minutes without moving the card.

10. The test is to be invalid after more than 25 minutes after dripping the analyte (sample) in to detection card.

[INTERFERENCE FACTORS]

No interference with the following substances has been detected: Benzocaine 150 mg/dL, Flixonase 50 MCG, Blood 5%, Tamiflu (Oseltamivir phosphate) 500mg/dL, Mucin 5 mg/mL, Pulmicort 0,5MG/m, NasoGel spray (Neil Med), Nutricius Biotin 300 mg, Neosynephrine-pois 100mg/mL, Methanol 150mg/dL, Afrin 0,5 mg/mL, Salbutamol 0,15 mg/mL, Robitussin expectorans 100 mg/mL, Mupirocin 0,15mg/dL, Vasoprin 100mg/tbl, Berodual N 10 ml, Psilo-Balsam 10mg/Gel, Alkalol Nasal Wash 10%, Dexromethorphan nutra essential 10mg/mL, Aniginal Nasal Mist (NeilMed), Fortecortin 4mg/tbl, Anoprin 100mg/tbl, Ventolin Inhaler N 100 mg, Symbicort Turbulizer 160mg/4,5mg.

[REPEATABILITY]

The repeatability and reproducibility study of three lots products from different laboratories and different dates, 10 Negative controls, 3 LOD controls and 1 CV control were used to study. The coincidence rate of intra-assay and inter-assay repeatability was 100%. The coincidence rate of intra-day and inter-day repeatability was 100%.

[LIMIT OF DETECTION - LOD]

The LOD for the Novel Coronavirus LOMINA SARS-CoV-2+Influenza A/B Antigen Test has been established using dilutions of a standard solution. The original concentration supplied is a concentration of SARS-CoV-2: 3.0 x 10³ TCID₅₀/mL; Influenza A: 1.0 x 10³ TCID₅₀/mL; and Influenza B: 1.0 x 10³ TCID₅₀/mL. A concentration was chosen between the last dilution to give 3 positive results and the first to give 3 negative results. Using the concentration which last dilution to give 3 positive results, the LOD was further refined with a 2-fold dilution series 10x. The last dilution demonstrating 100% positivity was then tested in an additional 20 replicates tested in the same way.

[SENSITIVITY]

Sensitivity: 98,71% (95% CI: 96,74%-99,65%)

Specificity: >99,9% (95% CI: 99,29%-99,98%)

Accuracy: 99,52% (95% CI: 98,77%-99,87%)

The results to SARS-CoV-2 virus of LOMINA SARS-CoV-2+Influenza A/B Antigen Test with correlation to Ct value of the positive samples were as follows:

CT Value	Diagnostic sensitivity LOMINA SARS-CoV-2+Influenza A/B Antigen Test	95% CI
PCR positive	98,70%	96,7-99,9%
PCR - Strongly positive (CT < 20)	100%	99,2-100%
PCR - Positive (CT 20-30)	99,2%	97,1-99,9%
PCR - Weakly positive (CT > 30)	84,6%	54,6-98,1%

[TEST RESULTS TO SARS-CoV-2 virus]

LOMINA SARS-CoV-2+Influenza A/B Antigen Test

PCR POSITIVE 307

PCR NEGATIVE 4

520

Sensitivity: 98,71% (95% CI: 96,74%-99,65%)

Specificity: >99,9% (95% CI: 99,29%-99,98%)

Accuracy: 99,52% (95% CI: 98,77%-99,87%)

The results to SARS-CoV-2 virus of LOMINA SARS-CoV-2+Influenza A/B Antigen Test with correlation to Ct value of the positive samples were as follows:

CT Value	Diagnostic sensitivity LOMINA SARS-CoV-2+Influenza A/B Antigen Test	95% CI
PCR positive	98,9%	96,9-99,9%
PCR - Strongly positive (CT < 20)	100%	90,3-100%
PCR - Positive (CT 20-30)	99,6%	99,3-100%
PCR - Weakly positive (CT > 30)	81,8%	48,2-97,7%

[TEST RESULTS TO Influenza A virus]

LOMINA SARS-CoV-2+Influenza A/B Antigen Test

PCR POSITIVE 284

PCR NEGATIVE 3

520

Sensitivity: 98,95% (95% CI: 96,98%-99,78%)

Specificity: >99,9% (95% CI: 99,29%-100%)

Accuracy: 99,63% (95% CI: 98,92%-99,92%)

The results to Influenza A virus of LOMINA SARS-CoV-2+Influenza A/B Antigen Test with correlation to Ct value of the positive samples were as follows:

CT Value	Diagnostic sensitivity LOMINA SARS-CoV-2+Influenza A/B Antigen Test	95% CI
PCR positive	98,9%	95,5-99,9%
PCR - Strongly positive (CT < 20)	100%	90,3-100%
PCR - Positive (CT 20-30)	99,6%	99,3-100%
PCR - Weakly positive (CT > 30)	81,8%	48,2-97,7%

[TEST RESULTS TO Influenza B virus]

LOMINA SARS-CoV-2+Influenza A/B Antigen Test

PCR POSITIVE 189

PCR NEGATIVE 3

520

Sensitivity: 98,44% (95% CI: 95,50%-99,68%)

Specificity: >99,9% (95% CI: 99,29%-100%)

Accuracy: 99,49% (95% CI: 98,77%-99,91%)

The results to Influenza B virus of LOMINA SARS-CoV-2+Influenza A/B Antigen Test with correlation to Ct value of the positive samples were as follows:

CT Wert	Diagnostische Sensitivität LOMINA SARS-CoV-2+Influenza A/B Antigen Test	95% CI
PCR positiv	98,70%	96,7-99,9%
PCR - Stark positiv (CT < 20)	100%	90,3-100%
PCR - Positive (CT 20-30)	99,2%	97,1-99,9%
PCR - Schwach positiv (CT > 30)	84,6%	54,6-98,1%

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